



Ústav fyzikální chemie J. Heyrovského, ČAV,
Praha 8, Dolejškova 3

Scanning Probe Microscopy I

Basic Techniques

Pavel Janda

Odd. elektrochemických materiálů

<http://www.jh-inst.cas.cz/>
<http://www.jh-inst.cas.cz/~janda>
pavel.janda@jh-inst.cas.cz

Microscopic methods/resolutions


OPT: Optical Microscopy

SNOM: Scanning Near Field Optical Microscopy


SEM: Scanning Electron Microscopy

HRTEM: High Resolution Transmission Electron Microscopy

STM, AFM:
Scanning Tunneling Microscopy,
Atomic Force Microscopy



Scanning Probe Microscopy





SPMs by information carriers







**Electron - tunneling microscopy STM/ECSTM,
conductive AFM CAFM (PFTUNA)**

Ion - electrochemical microscopy EC STM/EC AFM

Force interactions - atomic force microscopy AFM/ECAFM

Tunneling Microscopy

Binning, Rohrer, IBM, 1981, Nobel Prize 1986



Approximation of tunneling current

$$I_T \sim V_B f_{mTS}(V_B) \exp [-2z\sqrt{(2m\Phi_{ST}/\hbar^2)}]$$

$\hbar = h/2\pi$, $f_{mTS}(V_B)$...reduced Planck const.

$I_T(V_B)$ given by e-structure tip x sample,
z...distance tip-sample ($\sim 10^{-9}$ nm), V_B up to $\pm 1-2$ V, $I_T \sim nA - pA$

Tunneling barrier (distance) spectroscopy



Low bias $V_B = \text{const.}$: $(dI_T/dZ)/I_T \sim (2\sqrt{2m_e})/\hbar \sqrt{(\Phi_S + \Phi_T)}$

Φ_{Sample} , Φ_{Tip} local work functions, I_T tunneling current,
Z tip-sample distance, m_e e-mass

Realization:

modulated axial distance (Z-axis) VVVVV, recording $dI_T/dZ \Rightarrow \Phi_{S,T}$

Considering $\Phi_{\text{Tip}} \approx \text{const.} \Rightarrow$ lateral variation of barrier height \sim local Φ_S



-
-
-

Tunneling voltage spectroscopy



$V_B <$ work function tip-sample (~ 10 mV)

$dI_T/dV_B \approx$ local surface density of states real or derived from internal band arrangement

Modulation VVVVV V_B , recording I_T - V_B curves

Output: $d(\log I_T)/d(\log V_B)$ vs V_B


Yields: map of surface states (in UHV)
shows filling of states, ad-atoms
free (dangling) bonds...




I_T - V_B plots at single cryst. Si (UHV):
STM Tip passing over defect

[B. Persson, A. Baratoff, Phys.Rev.Lett. 59, 339]

Tunneling microscopy – topography uniform surface e-density



Tunneling microscopy – topography nonuniform surface e-density => artefacts




EC STM





Electrochemical Scanning Tunneling Microscopy

PJ


EC STM: Detection of tunneling currents in EC experimental setting



Probe/tip for EC STM



EC STM: Imaging Self-Assembled Monolayers on immersed surface




Aq. 0.1M Na₂HPO₄, pH 6

SAM - molecular ad-layer: mercaptopropionic acid (MPA) on Au111,
 $U_B = 300 \text{ mV}$ / $I_T = 100 \text{ pA}$

•
•
•

Tip-assisted „nanoprint“

STM tip-deposited nanoparticles Cu $d \geq 8$ nm, $h < 1$ nm





Atomic Force Microscopy

and force-based techniques



AFM

Force Interactions Range/mode

Long - magnetic, coulombic - **noncontact**

Medium - van der Waals (dipol-dipol, induced - dipol-nonpolar) **semicontact**

Short - binding (attractive) interactions, repulsive (deformation) – „**contact**“


„Derived AFM techniques“ ex-situ / in fluids

Conductive AFM (CAFM) - conductivity meas., **Tunneling AFM (TUNA) (!)**


Electrochemical AFM (EC AFM) – AFM imaging during EC experiments
+ tip-assisted „electro-lithographic“ techniques

Scanning Electrochemical Microscopy (SECM) – imaging based on
charge transfer reactions/Faraday currents (feedback)




AFM fundamentals



AFM fundamentals – force curve



AFM Probe: Tip and spring (*cantilever*)




material of tip and cantilever: Si, Si_3N_4

Source: Bruker



Scanning in selected regions of force curve

AFM in repulsive forces: Contact Mode

Maximum resolution, but maximum interaction tip-sample





AFM imaging in contact mode




Lateral XY~20 000x


Zeolite crystals




AFM imaging in contact mode



Graphite - HOPG





Mica -
muscovite




Oriented PTFE
molecules on Si wafer

Contact mode AFM - repulsive forces – indentation, nano-lithography





Contact mode AFM – attractive forces:

Adhesion
Friction



Attractive Lateral Forces (friction) for material analysis (LFM)


Contact Mode Lateral Forces



PTFE film on glass
-topography
-map of lateral forces
(friction)

•
•
•


Contact mode AFM – Dynamic Force Microscopy/Spectroscopy (DFM/S)



Energy dissipation -
“probe defectoscopy”:
elasticity (Y_M), plasticity

Tip in contact,
oscillating cantilever
 $f_{\text{drive}} \Rightarrow A_{\text{drive}}/A_{\text{response}}$






AFM DFS material analysis - Zr/Nb/ZrO₂



Attractive axial forces for „chemical“ analysis

Stretching molecules between cantilever tip and sample

Reference:



•
•
•

Techniques combining AFM and e-conductivity

Contact mode:

Conductive Force Microscopy (CFM)


Non-contact mode (!):
Tunneling AFM (TUNA)

P3


Conductive Force Microscopy (CFM) Tunneling AFM (TUNA)

Spreading Resistance imaging


Constant Force mode.



Copyright © NT-MDT, 2002



Carotenoid embedded in 1-docosanethiol attached to Au.
Current measured between biased Pt-coated AFM cantilever
and Au substrate.
Maximum current (I_m) vs. contact force.
[J. Phys. Chem. B 103 4006-4010 (1999)]



-
-
-

PeakForce Tapping QNM

Complete force curve is collected at each point of scan

For

- 3D Topography (contact, height) at chosen Peak Force
- Quantitative Nanoscale Mechanical analysis:
Mapping surface nanomechanical properties


adhesion
stiffness/ Y_M
deformation
dissipation

extracted simultaneously from various regions of force curve.

Feedback set by Peak Force (contact)Tapping



AFM PFQNM – Peak Force tapping

Figure 2.5f Force curve, Force vs. distance



Periodic contact with surface contact 3D (height) images


Recording of Force Curve



Source: Bruker


AFM PFQNM –interpretation of force curve: Deformation

Figure 2.5j Deformation




AFM PFQNM –interpretation of force curve: Adhesion

Figure 2.5c Adhesion on a PS+LDPE blend




Source: Bruker


AFM PFQNM –interpretation of force curve: Stiffness (Y_M /DMT)




AFM PFQNM – interpretation of force curve Dissipation

Figure 2.5h Dissipation (shaded area) in a polystyrene (PS) and Low-density polyethylene (LPDE) blend





Composite PS/LDPE
($Y_{PS} \sim 3$ GPa/ $Y_{LDPE} \sim 0.3$ GPa)
**AFM PFQNM -
interpretation
of force curve**



Deformation

Dissipation

Adhesion

Stiffness (Y_M , DMT-approx.)









Figure 2.5b DMT Modulus map of a PS+LDPE blend








Source: Bruker

AFM in liquids



PFQNM on soft layers - interpretation of force curve




No adhesion
to clean Si-wafer
(water-immersed)

Source: Lopez Mora et al - study on LBL films

Lipid bilayer on Si-wafer immersed
in aqueous media – delay on rising part
of force curve: compression,
penetration resistance (PFQNM)


•
•
•

AFM imaging of “zero-mass” nano-objects: Gaseous nanobubbles on immersed surfaces




In-situ AFM *tapping*, amplitude
image of surface nanobubbles
on water-immersed HOPG.
Hydrophilized tip and cantilever

AFM-PFQNM of gaseous nanodomains - interpretation of force curve






Gaseous nanodomains on water-immersed surface – hydrophilized tip




Hana Tarabkova and Pavel Janda: Effect of Graphite Aging on Its Wetting Properties and Surface Blocking by Gaseous Nanodomains, Langmuir 2023, 39, 14154–14161

In-situ AFM imaging of electrochemical experiment




Cu deposit on HOPG basal plane removed by SPM tip

AFM in Scanning Electrochemical Microscopy (SECM)



AFM-SECM tip/electrode
insulated by electrophoretic
coating, tip $r \sim 130$ nm

Feedback: Farad. current




•

•


•

AFM in EC: ECAFM/lithography tip-assisted electrodeposition of metal in droplet




Au electro-deposition: conversion Au(III) (ink) => Au(0)
2-electrode system

Constant Force mode




ECAFM/lithography: tip-assisted local electrochemical deposition/dissolution of oxide





Field intensity 10^7 V/cm =>
water dissociates $H_2O \Rightarrow OH^- + H^+$

2-electrode system



Water meniscus on hydrophilic
and hydrophobic substrate







Semi(contact) mode - tapping

Acoustic/mgt. drive

Mechanical oscillator
oscillation in resonance
input parameters:
 $f_{rez} A_{sp}$ (~ 20 nm)

output parameters
 $A, \Delta f, \Delta\theta, d$ (deflect.)

"V"cantilever ($k = 0.58$ N/m)



$$md^2z/dt^2 = -kz + (m\omega_0/Q)dz/dt + F_{tip-sample} + F_{drive} \cos\omega t$$

Hook


Dissipation
(dumping)

Interaction
tip-sample


$$\omega_0 = \sqrt{k/m}$$

Semicontact mode (tapping) phase imaging

Enhances edges of (nano)structure boundaries.
Less affected by height fluctuation





Nanograins ZrO_2 imaged by AFM
tapping phase mode (top-view)










Identical location imaged by 3D topography
and amplitude mode (top-view)


AFM artefacts



AFM artefacts




AFM artefacts: Self-image of AFM tip




• • •

Special: Dual function of micropipette: AFM scanning tip and micro-dispenser




Nano-lithographic AFM with scanning micropipette




F. Iwata


Special: AFM *in vivo* - Scanning Force Endoscope



Diagnosing cartilage diseases at an early stage



M. Stolz et al., Biophys. J. 2004; 86 3269-3283



Institute of Microtechnology
University of Neuchâtel


R. Imer
T. Akiyama
N.F. de Rooij
U. Staufer

M.E. Müller Institute (MSB)
University of Basel

M. Stolz
U. Aebi



Ústav fyzikální chemie J. Heyrovského, ČAV,
Praha 8, Dolejškova 3
Odd. Elektrochemických materiálů



<http://www.jh-inst.cas.cz/~janda>
pavel.janda@jh-inst.cas.cz

PJ

JT

